

BD Unopette™ System Test

REF 365854
REF 365855

IVD

**Sistema
Système test
Systemtest
Test**

Pages

1-2	EN	English
3-4	ES	Español
5-6	FR	Français
7-8	DE	Deutsch
9-10	IT	Italiano

Symbols

EN Keep Away From Sunlight
ES No exponer a la luz solar
FR Protéger de la lumière solaire
DE Vor Sonnenlicht schützen
IT Non esporre alla luce solare

EN Manufacturer
ES Fabricante
FR Fabricant
DE Hersteller
IT Produttore

EN Authorized Representative
ES Representante autorizado
FR Représentant autorisé
DE Bevollmächtigter Vertreter
IT Rappresentante autorizzato

EN Fragile, Handle With Care
ES Frágil, manipular con precaución
FR Fragile, manipuler avec précautions
DE Zerbrechlich, mit Sorgfalt behandeln
IT Fragile, maneggiare con cura

EN This End Up
ES Este lado hacia arriba
FR Se terminez vers le haut
DE Diese Seite ist oben
IT Estremità superiore

LOT Batch Code

REF Catalog Number

Use By

Do Not Reuse

Consult Instructions For Use

Caution, Consult Accompanying Documents

Upper Limit of Temperature

In Vitro Diagnostic Medical Device

Product Information

Intended Use

BD Unopette™ System Test is a stable *in vitro* diagnostic reagent system for the enumeration of leukocytes and platelets in whole blood.

Summary

This method is based on the procedure developed by Brecher and Cronkite.¹ The ammonium oxalate diluent is modified by the addition of Sorensen's phosphate buffer to maintain the pH of the reagent, and thimerosal to act as an antibacterial agent.

Principle

Whole blood is added to the diluent, the tonicity of which lyses red cells but preserves platelets, leukocytes and reticulocytes. When erythrocytes are completely lysed, the solution will be clear red and counting can proceed. The diluted blood is placed in a hemacytometer according to accepted technique. Cells are allowed to settle for ten (10) minutes before leukocytes and platelets are counted.

Under 100X magnification using bright-light microscopy, leukocytes appear refractile. At this magnification some red cell stroma will be evident but will not interfere with the counting procedure.

Under 430X magnification using bright-light microscopy, platelets appear oval or round and frequently have one or more dendritic processes. Under the same magnification using phase microscopy, they appear morphologically the same as with bright-light microscopy but are uniformly opaque with a pink or purple sheen.

Leukocyte counts obtained with BD Unopette™ System Test compare favorably with those obtained using the Thoma White Cell Pipette with 3% acetic acid. Platelet counts obtained with BD Unopette™ System Test compare favorably with those obtained using the Thoma Red Cell Pipette with Rees-Ecker diluent. Platelet counts obtained with BD Unopette™ System Test also compare favorably with results obtained with the Coulter Counter®* Model F electric counting method.

BD Unopette™ System Test

1. BD Unopette™ reservoir containing 1.98 ml of 1% Buffered Ammonium Oxalate
2. BD Unopette™ capillary pipette 20 µl capacity

Warnings:

Moisture is added to BD Unopette™ reservoir jars to maintain shelf life. Thimerosal, a preservative added to both the jar and the reservoir reagent, may cause skin irritation. If reservoir is squeezed too hard, the specimen may be expelled through the top of the overflow chamber, resulting in contamination of the exposed skin. All specimen handling, diluting, and handling of the filled hemacytometer should be performed wearing personal protective equipment (i.e. gloves, lab coat, eye protection, etc).

Dilution Ratio

Sample to total volume 1:100

Reagent Stability

The reagent system is stable if stored below 30°C (86°F), protected from sunlight, and used before the expiration date under the stated conditions of the test. Do not use if reagent solution is not clear.

Equipment

1. Microscope equipped for 100X and 430X magnification. Phase contrast for platelet counts is optional.
2. Hemacytometer. Neubauer hemacytometer recommended.
3. Petri dish. Bottom should be lined with moist filter paper.

Specimen

Free-flowing capillary or thoroughly mixed anticoagulated venous blood. EDTA is the anticoagulant of choice. Use of an appropriate BD Vacutainer® Tube with EDTA is recommended for collection of venous specimens. (See Technical Notes 1 and 2). The diluted specimen is stable for three (3) hours at room temperature.

Interferences

Bacteria and other debris can interfere with platelet counts. However, bacteria, dirt and crystals can usually be properly identified by their refractility and lack of pink or purple sheen under phase microscopy. Care must be taken to differentiate platelets from extrinsic particles under bright-light microscopy.

Storage Instructions

If count cannot be performed immediately after blood is diluted, store reservoir at room temperature. Perform count within three (3) hours of making dilution.

*Coulter Counter® is a registered trademark of Beckman Coulter, Inc. Fullerton, California.

Procedure

1. Puncture Diaphragm

Using the protective shield on the capillary pipette, puncture the diaphragm of the reservoir as follows:

- a. Place reservoir on a flat surface grasping reservoir in one hand, take pipette assembly in other hand and push tip of pipette shield firmly through diaphragm in neck of reservoir, then remove.

- b. Remove shield from pipette assembly with a twist.

2. Add Sample

Fill capillary with whole blood and transfer to reservoir as follows:

- a. Holding pipette almost horizontally, touch tip of pipette to blood. See adjacent diagrams. Pipette will fill by capillary action. Filling is complete and will stop automatically when blood reaches end of capillary bore in neck of pipette.

- b. Wipe excess blood from outside of capillary pipette, making certain that no sample is removed from capillary bore.

- d. Cover opening of overflow chamber of pipette with index finger and seat pipette securely in reservoir neck.

- c. Squeeze reservoir slightly to force out some air. Do not expel any liquid. Maintain pressure on reservoir.
- e. Release pressure on reservoir. Then remove finger from pipette opening. Negative pressure will draw blood into diluent.

- f. Squeeze reservoir gently two or three times to rinse capillary bore, forcing diluent up into, but not out of, overflow chamber, releasing pressure each time to return mixture to reservoir. See Caution Section.
- g. Place index finger over upper opening and gently invert several times to thoroughly mix blood with diluent.
- h. Let stand for ten (10) minutes to allow red cells to hemolyze. Leukocyte and platelet counts should then be performed within three (3) hours.

3. Charge Hemacytometer

Mix diluted blood thoroughly by inverting reservoir (see 2g above) to resuspend cells.

- a. Convert to dropper assembly by withdrawing pipette from reservoir and reseal securely in reverse

- b. To clean capillary bore, invert reservoir, gently squeeze sides and discard first three or four drops.

- c. Carefully charge hemacytometer with diluted blood by gently squeezing sides of reservoir to expel contents until chamber is properly filled.

- d. Place hemacytometer on moistened filter paper in Petri dish. Cover Petri dish and allow to stand ten (10) minutes to permit cells to settle. (Moistened filter paper retards evaporation of diluted specimen while standing.)

4. Count and Calculate

A leukocyte count is performed with a Neubauer hemacytometer as follows:

- a. Under 100X magnification, leukocytes are counted in all nine (9) large squares of the counting chamber.
- b. Add 10 percent of count to total number of cells counted. (See Technical Note 4)
- c. Multiply this figure by 100 to get total leukocyte count. Example: If 70 cells are counted, total count is: (70 + 7) X 100 = 7700 leukocytes/cu mm.

A platelet count is performed with a Neubauer hemacytometer as follows:

- a. Under 430X magnification using bright-light or phase microscopy, platelets are counted in all 25 small squares within the large center square.
- b. Multiply number of platelets counted by 1000 to get total platelet count.

Example: If 150 platelets are counted, total count is: 150 X 1000 = 150,000 platelets/cu mm.

Note: If another hemacytometer is used, calculations appropriate to the particular type should be employed.

Limitation of Procedure

A highly elevated leukocyte or platelet count may make accurate counting difficult. In either instance, a secondary dilution should be made. When calculating the total count, adjust the formulas to allow for secondary dilution.

Expected Performance

1. Comparison of platelet counts of 526 specimens of blood, performed in duplicate by Coulter Counter® Model F and BD Unopette™ System Test.

	Platelets per cu mm	
	BD Unopette™ System Test	Coulter Counter® Model F
Mean	244,900	207,700
Precision to one standard deviation by comparison of duplicates	± 15,400	± 15,700
Coefficient of variation	6.3%	7.6%

2. Precision of leukocyte counts and platelet counts using BD Unopette™ System Test.

Leukocyte Count

Determination (n)	40
Mean	5221 Leukocytes / cu mm
Standard deviation	± 382 Leukocytes / cu mm
Coefficient of variation	7.3%

Platelet Count

Determination (n)	18
Mean	303,000 Platelets / cu mm
Standard deviation	± 18,000 Platelets / cu mm
Coefficient of variation	5.9%

Technical Notes

1. Platelets become activated by surface contact and will aggregate. If samples are collected into evacuated tubes, tubes must be inverted 8-10 times immediately after filling. Delay in mixing may result in platelet clumping and erroneous results.
2. Excessive tissue trauma during collection of blood may also result in platelet clumping and reduce the accuracy of the platelet count. When clumps of platelets are observed in a hemacytometer, another sample of blood must be collected.
3. Cells and diluent must be adequately mixed and counting chambers should be properly filled if errors in manual counting procedures are to be avoided.
4. This step simplifies the calculation which actually entails dividing the number of cells by the number of squares counted and multiplying by 10 to correct for the depth of the chamber.

References

1. Brecher G, Cronkite EP: Morphology and enumeration of human blood platelets, J Appl Physiol 3:365 (Dec) 1950.

Caution:

Handle all biologic samples and blood collection products according to the policies and procedures of your facility.

All glass has the potential for breakage; take precautionary measures during handling.

Obtain appropriate medical attention in the event of any exposure to biologic samples (e.g., through a puncture injury) since samples may transmit viral hepatitis, HIV (AIDS), or other infectious diseases.

Discard all products contaminated with blood in biohazard containers approved for their disposal.

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