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





### Sets With Multiple Sample Luer Adapter

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

### Sets With Female Luer Luer Adapterなし





### Withdrawing Safety-Lok™ Blood Collection Sets

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# VACUTAINER® Brand Safety-Lok™ Blood Collection Sets

NOTE: Venipuncture should be performed wearing gloves.

1. Peel apart package and remove set.

2. **WITH MULTIPLE SAMPLE LUER ADAPTER:** Thread the Luer Adapter into Holder (Figure 1). A Becton Dickinson VACUTAINER® Brand Needle Holder or a Becton Dickinson Adapter of an appropriate size must be used with this product to ensure proper function. Check to ensure that the Female Luer (A) is securely attached on the Male Luer (B).

**WITH FEMALE LUER:** Remove tip cap (Figure 1A). Assemble syringe (Figure 2A).

**Note:** Filling the tubes from a hypodermic syringe while the stopper is in place is not recommended. Forcefully depressing the syringe plunger without removing the stopper can create positive pressure in the tube causing the stopper and specimen to fly out with explosive force.

3. Utilizing tourniquet, if appropriate, identify site for venipuncture.

4. Prepare venipuncture site with appropriate antiseptic. **Allow site to dry.**

5. Remove needle sheath (Figure 2 or 2A) and perform venipuncture by holding wings as shown (Figure 3 or 3A). Do not grasp the yellow safety shield to insert device.

6. **WITH MULTIPLE SAMPLE LUER ADAPTER:** Collect into an evacuated tube by placing tube into holder, then puncture the stopper by pushing the tube and replace with additional tubes as appropriate. (See order of draw). Do not leave tourniquet in place for more than 2 minutes. Remove as soon as blood appears in the first tube. A draw loss of approximately 0.5 ml will result when collecting the initial tube. **Tubes containing additives should be mixed as soon as they are filled.**

**WITH FEMALE LUER:** Use Safety-Lok™ Blood Collection Set according to the Standard Operating Procedures of your institution. If syringe is used for blood collection, see instructions for use of Blood Collection Tubes and NCCLS Approved Standard H3-A3, Vol. 11, No. 10, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.

7. **Withdrawing Safety-Lok™ Blood Collection Sets:** When final tube is filled, remove the Safety-Lok™ Blood Collection set by:

- Grasping the translucent yellow safety shield grip area with the thumb and index finger while at the same time grasping the tubing (Figure 4), or
- Grasp either one wing or both wings and withdraw (Figure 4A). Apply digital pressure to site using sterile gauze pad.

8. **Safety-Lok™ Shield Activation:**

- **One-Hand Technique.** Hold tubing in hand and advance translucent yellow safety shield with thumb and index finger (Figure 5) until a click is heard indicating that the needle is completely retracted and safety shield is locked in place over needle tip (Figure 8), or...
- **Modified One-Hand Technique.** Apply pressure to site using your fingers as shown (Figure 4 or 4A). Withdraw blood collection set by grasping the translucent yellow safety shield grip area with the thumb and index finger. With opposite hand, grasp tubing between thumb and index finger while pushing the yellow safety shield forward (Figure 6) until a click is heard indicating that the needle is completely retracted and the safety shield is locked in place (Figure 8), or...

- **Two-Hand Technique (Figure 7).** Grasp either wing with one hand (A) and grip area of the yellow safety shield base with the other hand (B). Slide the wings back into the rear slot of the safety shield, until a click is heard indicating that the needle is completely retracted and locked into place over needle tip (Figure 8).

**Do not grasp safety shield by body; this will interfere with safety shield activation.**

9. Label all specimens.

10. **Safety-Lok™ Blood Collection Set Disposal:** Dispose of the holder and Blood Collection Set into the nearest sharps container according to state and local regulations. If separating holder from Blood Collection Set, use caution as the Multiple Sample Luer Adapter is unprotected. Dispose of all used materials in appropriate container and wash hands.

#### RECOMMENDED ORDER OF DRAW:

1. Tubes for sterile samples.
2. Tubes without additives.
3. Tubes for coagulation studies (e.g., citrate).
4. Tubes with other additives (e.g., heparin, EDTA).

**SST® Tubes, VACUTAINER® Clot Activator Tubes, and VACUTAINER® PLUS Tubes are considered additive tubes.**

#### CAUTION:

Handle all biologic samples and blood collection "sharps" (lancets, needles, and blood collection sets) in accordance with the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since the samples may transmit HBV (Hepatitis), HIV (AIDS), or other infectious diseases.

Utilize any built-in used needle protector, if the blood collection device provides one.

Discard all blood collection "sharps" in biohazard containers approved for their disposal.

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