ARCTIC SUN™

TEMPERATURE MANAGEMENT SYSTEM

CALIBRATION TEST UNIT (CTU) OPERATOR'S MANUAL

TABLE OF CONTENTS

Section 1 — Introduction	4
1.1 Use of the Calibration Test Unit (CTU) Operator's Manual	4
1.2 System Description	4
1.3 Symbols and Standards	4
1.4 Environmental Conditions	5
1.5 General Warnings	5
Section 2 — Instructions for Use	5
2.1 Overview	5
2.2 Initiating a Calibration or Test	5
Section 3 — Theory of Operation	6
3.1 Introduction	6
3.2 Basic CTU Operation	6
3.3 Controls	7
3.4 Connections	7
3.5 Understanding the Test and Calibration Sequence	7
3.5.1 Calibration Factor Longevity	7
3.5.2 Initiation	7
3.5.3 Pre-Warm and Flow Check	7
3.5.4 Patient Temperature Channel Checks	7
3.5.5 Temperature Out Check	8
3.5.6 Water Temperature Check	8
3.5.7 Temperature Out Calibration	8
3.5.8 Check Calibration	9
Section 4 — Maintenance	9
4.1 Recommended Maintenance	9
4.2 Calibration	9
4.3 Technical Assistance/Service	9
4.4 Disposal	9
APPENDIX A: Specifications	10
APPENDIX B: Error Codes	11

Section 1 – Introduction

1.1 Use of the Calibration Test Unit (CTU) Operator’s Manual

The Calibration Test Unit provides an effective and easy way of checking and calibrating, if necessary, the ARCTIC SUN™ Temperature Management System. This Manual provides a detailed description of the CTU, its use, construction, routine maintenance and calibration. It is important that the CTU be properly maintained and serviced in accordance with the instructions described in this manual. Only trained personnel should use and service this product.

The CTU is intended to be used in association with the ARCTIC SUN™ Temperature Management System Service Manual. The Service Manual provides detailed information about the ARCTIC SUN™ Temperature Management System and how the CTU can be used to troubleshoot problems in addition to its normal test and calibration functions.



1.2 System Description

The CTU is, for the most part, a passive device that provides stimulus and feedback to the ARCTIC SUN™ Temperature Management System to enable it to test and calibrate itself. The sole active function of the CTU is to provide an external heat source to raise the circulating water temperature for troubleshooting under the guidance of Medivance technical support personnel.

1.3 Symbols and Standards

The CTU bears the following symbols affixed to the system:

Table 1.1. CTU Symbols

	For the safe and effective use of this device, the operator must consult the accompanying documents prior to use.
	Per ETL Intertek, models of the ARCTIC SUN™ Temperature Management System CTU that bear the ETL Monogram have been certified to CSA STD C22.2 No. 61010-1 and conform to UL STD 61010-1.
	Indicates high temperature part or component. The maximum temperature of this internal component allowed by the protective system is also listed.
	Indicates electrical hazard.
	Indicates Earth Ground.
	Identifies the heater fuse.
	Identifies the storage temperature range.
	Identifies the storage relative humidity range.
	Indicates the date of manufacture.
	Denotes calibration data.
	Denotes a push-button, momentary switch.
	Manufacturer
	Catalog number
	Indicates the Authorized Representative in the European Community.
	The CTU must be disposed of properly. DO NOT dispose of the unit into the garbage.
	Authorized Representative in Switzerland
	Packaging unit
	Serial number

1.4 Environmental Conditions

- Ambient Temperature Range:
 - o Operating Temperatures: 18°C to 24°C (65°F to 75°F)
 - o Transport and Storage Temperatures: -30°C to 50°C (-20°F to 120°F)
- Ambient Humidity Range:
 - o Operating Humidity: 5% to 95% relative humidity, non-condensing
 - o Transport and Storage Humidity: 5% to 95% relative humidity, non-condensing
- Ambient Atmospheric Pressure Range:
 - o Operating: 375 mmHg to 795 mmHg
 - o Transport and Storage: 375 mmHg to 795 mmHg

NOTE: If the Calibration Test Unit is used at higher room temperatures, the test system’s cooling capability and accuracy may be degraded.

- Ingress Protection Rating IPX0 for ordinary usage.

1.5 General Warnings

- This equipment is to be used only as described in this manual. Failure to do so may result in damage to the device.
- The CTU must not be used on equipment while connected to a patient.
- Do not use the Calibration Test Unit in the presence of flammable agents because an explosion or fire may result.
- Anyone performing the procedures must be appropriately qualified.
- Carefully observe the Calibration Test unit for fluid leaks before and during use. Leakage may result in lower flow rates or an electrical shock hazard. Discontinue use immediately if leakage is observed.
- For indoor use only.

Section 2 – Instructions for Use

2.1 Overview

The entire test and calibration process takes approximately 90 minutes. Operator assistance is required only in certain stages:

- Stage 1: Warm-up, Flow Check, and Inlet Pressure Calibration: 12 minutes, automatic.
- Stage 2: Patient Temperature Calibration: 5 minutes, operator assistance required.
- Stage 3: Water Temperature Calibration: 18 minutes, operator assistance required (at completion).
- Stage 4: Water Temperature Calibration, Heater Check: 33 minutes, operator assistance required (at completion).
- Stage 5: Temperature Out Calibration: 25 minutes, operator assistance required (at completion).

2.2 Initiating a Calibration or Test

- A. Replace the fluid delivery line with CTU.
- B. Connect the Blue circular connector labeled “PT1” to the connector Patient Temperature 1 (larger thermometer and patient symbol).
- C. Connect the Blue circular connector labeled “PT2” to the connector Patient Temperature 2 (smaller thermometer and patient symbol).
- D. Connect the Black circular connector labeled “TO” to the connector labeled “TEMP OUT”.
- E. Power on the ARCTIC SUN™ Temperature Management System Control Module.
- F. Press the Advanced Setup button on the Therapy Selection screen.
- G. Press the Start button next to Calibration on the Advanced Setup screen.
- H. Select Sensor Calibration to perform a calibration or calibration check of the ARCTIC SUN™ Temperature Management System.

- I. Adjust the values A to H on the CTU calibration screen to match the values on the CTU label. Press Continue when complete.
- J. Select either New Calibration or Check Calibration and follow the on-screen instructions.


Figure 2.1. ARCTIC SUN™ Temperature Management System Control Module Calibration Screen


Figure 2.2. ARCTIC SUN™ Temperature Management System Control Module Calibration Screen


Figure 2.3. ARCTIC SUN™ Temperature Management System Control Module Calibration Screen

Figure 2.4. ARCTIC SUN™ Temperature Management System Control Module with fluid delivery line removed

Figure 2.6. CTU switch 1


Figure 2.7. CTU switches 2 and 3

Figure 2.5. ARCTIC SUN™ Temperature Management System Control Module with CTU attached

Figure 2.8. CTU calibration data

Section 3 – Theory of Operation

3.1 Introduction

The Calibration Test Unit (CTU) is an instrument that is used to calibrate and verify operation of the ARCTIC SUN™ Temperature Management System Control Module. It is a hand sized unit that connects directly to the fluid delivery line port, taking the place of the fluid delivery lines and pads. The CTU can be used to semi-automatically test and calibrate the Control Module's two patient temperature measurement inputs, the temperature out function, and the internal water temperature measurement channels. It can also be used to automatically verify the Control Module's flow rate and inlet pressure measurements.

3.2 Basic CTU Operation

The CTU is used as an external stimulus and measuring device by software embedded in the Control Module. By testing the subsystems of the CTU in sequence, operation of the Control Module can be verified and calibrated.

The CTU contains precision resistors that simulate several patient temperatures. After a warm-up period, the simulated temperatures are first applied to the system patient temperature inputs to temporarily calibrate each channel. This calibration is temporary in the sense that it is used to perform the balance of the testing. It is not saved by the Control Module until the calibration process is complete.

Once the patient temperature channel's accuracy has been temporarily established, the CTU's temperature probe is connected to the Control Module's Patient Temperature 1 input. The Control Module then circulates water through the CTU at various temperatures. The patient temperature channel can then be used to check and calibrate the accuracy of the internal water temperature sensors and circuits.

The fluid path in the CTU is designed to allow a specific flow, which was selected to provide this flow rate at a given pressure. Since the ARCTIC SUN™ Temperature Management System controls the inlet pressure at a constant value, a crosscheck of flow and pressure measurement accuracy can be obtained. The flow rate and inlet pressure offset are checked, and the inlet pressure offset is zeroed during the initial warm-up period.

The CTU contains heaters that can be used to raise the system water temperature for troubleshooting. The test unit heaters are thermally and electrically protected, and are only powered while pressing the Heater Power Switch.

3.3 Controls

A rotary switch (S1), two slide switches (S2 and S3), and a heater switch (S4) comprise the CTU controls. The Output Select Switch, S1, is used to select the temperature simulator to be applied to the patient temperature channel. Six switch positions represent specific temperatures as follows:

- Position A is "open" or "no connection"
- Position B is approximately 10.8°C
- Position C is approximately 25.5°C
- Position D is approximately 34.5°C
- Position E is approximately 40.0°C
- Position F is a "shorted" connection.

The precise values of the simulators listed above are determined by a factory calibration process and are listed on a calibration sticker on the CTU.

The slide switch, S2, is used as follows:

- Position A connects Patient Temperature Channel 1 to the rotary switch and Patient Temperature Channel 2 to the Temperature Out Channel.
- Position B connects Patient Temperature Channel 1 to the CTU temperature sensor and Patient Temperature Channel 2 to the rotary switch.
- Position C connects Patient Temperature Channel 1 to the Temperature Out Channel.

The slide switch, S3, is used as follows:

- Position A connects the internal thermistor in the water flow path to the patient temperature circuit.
- Position B connects a parallel resistance with the internal thermistor, shifting the temperature up at the low water temperature setting, allowing it to be read by the patient temperature circuit.

The heater power switch, S4, is a momentary, normally open switch that applies AC power to the two 250W heaters when pressed and held in the ON position. The AC power cable must be connected to the CTU and an appropriate power receptacle to power the heaters.

3.4 Connections

A hydraulic connection and three electrical connections must be made to the control module so that the unit will function correctly in conjunction with the control module's software functions.

- The hydraulic connection is a custom fitting which mates directly to the control module. The seal must be tight with no leaks to ensure a proper test.
- The three electrical connections are labeled 'PT1', 'PT2', and 'TO'.
 - o PT1 corresponds to the primary patient temperature receptacle on the control module.
 - o PT2 corresponds to the secondary patient temperature connector.
 - o TO corresponds to the temperature out connector on the control module.

- TO has a different keying scheme from PT1 and PT2, and therefore cannot be interchanged. However, care must be taken when connecting PT1 and PT2 because they can be interchanged. If this does occur, the control module will indicate an error on the user interface display.
- The PT1, PT2, and TO electrical cables do not supply a voltage or current. They are passive connections with the voltage or current being supplied by the control module.

3.5 Understanding the Test and Calibration Sequence

The Calibration Test Unit and the ARCTIC SUN™ Temperature Management System Control Module work together to perform the tests and calibration in a semi-automated sequence. This section gives a step by step description of how this is accomplished. The calibration process may be discontinued at any time by pressing the Cancel button. Power must be cycled after a New Calibration process or Check Calibration process has begun to resume normal operation.

3.5.1 Calibration Factor Longevity

In order to perform this process accurately, various parameters of the Control Module must first be calibrated to be used in subsequent measurements. It should be noted that this calibration is temporary and it applies only during the Test and Calibration Sequence. This data is saved only at completion of the New Calibration sequence.

3.5.2 Initiation

The Control Module will enter Test and Calibration Mode only if all of the following are true:

- Patient Temperature 1 is between 10°C and 12°C.
- Patient Temperature 2 is less than 10°C.
- The reservoir level is 4 or greater.

These unique conditions are required to ensure that the system is connected to the CTU, not a patient. If one of the criteria above is not met, the user is prompted to rectify it.

3.5.3 Pre-Warm and Flow Check

The system pre-warm and flow check begins when the Test and Calibration Mode is initiated. This process starts with an automatic 10-minute pre-warm period. During this time, the unit runs in Manual Mode with a 28°C water target temperature and a countdown timer displays how many minutes remain in the pre-warm step. Approximately 55 seconds after initiating pre-warm, a bypass mode flow check is also performed.

For successful completion of pre-warm to occur, the following criteria must be met after 10 minutes:

- The bypass mode flow check flow rate must be greater than the non-bypass flow rate + 500ml/min.
- The average flow rate must be within ± 600 ml/min. of the CTU H value.
- The inlet pressure average must be within 1.0 psi of -7.0 psi.
- The bypass valve must be closed.

Once these conditions have been met, the system performs a purge for approximately one minute to prepare for checking the inlet pressure offset. The system performs an internal test of the patient temperature channels when purging is complete. If successful completion of the pre-warm / flow check does not occur within 12 minutes of its initiation, a time-out error is indicated.

At this point, the inlet pressure offset is calculated after the pressure stabilizes within 0.01 PSI. If the offset is greater than ± 0.5 PSI, an error message is displayed.

3.5.4 Patient Temperature Channel Checks

Operator assistance is needed to perform the Patient Temperature 1 (PT1) checks. The operator is prompted to change CTU Switch 1

settings which sequences the system through measurement of a set of precision resistors that simulate a patient temperature sensor at various temperatures. Each check is completed when the following criteria is met:

- The difference between patient temperature average and the respective CTU temperature value is within 0.8°C.
- The patient temperature is stable within 0.01°C.

If any check cannot meet the above criteria within 60 seconds, an error message is displayed. The check may be retried by selecting Continue.

Temporary calibration parameters are calculated based on the data obtained from Check 1 (CTU value A) and Check 4 (CTU value D). These parameters are then applied to the average values obtained for all four checks. These values are then checked against the respective CTU values that had been entered. The measured values must be within 0.14°C of the CTU values for successful completion. An error message is displayed if the criteria is not met.

The process is identical for Patient Temperature 2 (PT2) with the exception of the S2 position.

3.5.5 Temperature Out Check

The Temperature Out Check is only performed when Check Calibration has been selected. The Temperature Out Check is very similar to the patient temperature checks in that the same four sensor simulators are used. The operator is prompted to sequence CTU Switch 1 through a PT1 measurement of each, while PT2 measures the Temperature Out. (Temperature Out is an “echo” of Patient Temperature 1.) Obviously, this test requires precise measurement by the patient temperature channels. Therefore, the patient temperature checks must first be successfully completed.

Each check is completed when the following criteria is met:

- The difference between the Temperature Out average (as measured by PT2) and the respective CTU temperature value is within 1.0°C.
- The difference between the PT1 average and Temperature Out average (as measured by PT2) is within 0.8°C.
- The Temperature Out (as measured by PT2) is stable within 0.01°C

If any check cannot meet the those criteria within 75 seconds, an error message is displayed.

3.5.6 Water Temperature Check

The system prompts the user to change the S3 position at the beginning of the low temperature setpoint and then back after the low temperature setpoint has completed. The system runs in Manual Mode for the duration of the test. As the system runs, the water temperature is controlled at various setpoints while PT1 measures the temperature of the water flowing through the CTU. This reference temperature is then used to check the system’s internal water temperature sensors. The heater elements are checked, as well as, heating and cooling capacity tests.

Four segments which normally take about 15 minutes each comprise the Water Temperature Check. Check 1 takes place at approximately 6°C, Check 2 at approximately 28°C, Check 3 at approximately 40.5°C, and Check 4 at approximately 38°C. During this process, Patient Temperature 1, Outlet Water Temperature 1 (T1), Outlet Water Temperature 2(T2), Inlet Water Temperature (T3), Chiller Temperature (T4), the inlet pressure, and the flow are continuously averaged and monitored. Each check point requires the following to complete successfully:

- The average T2 must be within 2.0°C of the corrected reference temperature (i.e. the CTU probe as measured by PT1).
- The corrected reference temperature must momentarily cross over the desired cal. temperature (6°C, 28°C, and 40.5°C) and then cross back over again. The only exception is check 4 where the reference temperature must only drop 6 degrees below 38°C.
- The Bypass Valve must be closed.

If a check cannot be completed in 30 minutes, a message is displayed.

Again, this test requires precise measurement by the patient temperature channels. Therefore, the patient temperature checks must first be successfully completed. Further, in the measurement of water temperature, flow is also important. Therefore, the inlet pressure check must also have been completed as suction is the driving force in creating flow.

Newly calibrated temporary values are checked against the corrected average probe temperature values obtained during the checks. Both T1 and T2 readings must be within ± 0.54°C of the probe at 6°C and 28°C and within ± 0.34°C of the probe at 40°C, the T3 sensor must be within ± 0.54°C of the probe reading at 6°C, 28°C, and 40°C, and the T4 sensor must be within ± 0.54°C (checked at the 6°C temperature only). An error message is displayed if the criteria is not met.

3.5.7 Temperature Out Calibration

The final stage of calibration uses PT1 to measure the Temperature Out as it is stepped through the range of possible values. This process takes about 28 minutes. Upon completion, the calibration values are saved.

The information saved during the save process includes;

- New calibration factors generated during the various checks.
- CTU values (A through H) entered at the beginning of the test.
- Calibration Hours value which was set to zero prior to the save.

The stored CTU values are set to defaults prior to conducting the first calibration with the CTU. Additionally, stored CTU values will automatically be set to defaults if the validity check value for any of the CTU values is incorrect when retrieved from memory.

The default CTU values are:

- A. 10.87°C
- B. 25.54°C
- C. 34.52°C
- D. 39.98°C
- E. 0.00°C
- F. 0.00°C
- G. 0.00°C
- H. 2900

3.5.8 Check Calibration

The final stage of the Check Calibration process displays the values resulting from the check along with a pass or fail indication. The pass / fail criteria for the parameters shown on the result screen are indicated in Table 3.1.

Test	Criteria
Flow Rate	$300\text{ml/min} \leq \text{Flow} \Delta \leq 300\text{ml/min}$
Bypass Flow Rate	$\text{Bypass Flow} \Delta \geq 500\text{ml/min}$
Inlet water pressure zero offset	$-0.3 \text{ psi} \leq P \leq 0.3 \text{ psi}$
Patient Temperature 1 (Control 1) PT1 (4 checks)	$-0.3^{\circ}\text{C} \leq 10^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 25^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$ $-0.1^{\circ}\text{C} \leq 35^{\circ}\text{C} \Delta \leq 0.1^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$
Patient Temperature 2 (Monitor 2) PT2 (4 checks)	$-0.3^{\circ}\text{C} \leq 10^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 25^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$ $-0.1^{\circ}\text{C} \leq 35^{\circ}\text{C} \Delta \leq 0.1^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$
Patient Temperature Echo Out (Out) (4 checks)	$-0.2^{\circ}\text{C} \leq 10^{\circ}\text{C} \Delta \leq 0.2^{\circ}\text{C}$ $-0.2^{\circ}\text{C} \leq 25^{\circ}\text{C} \Delta \leq 0.2^{\circ}\text{C}$ $-0.2^{\circ}\text{C} \leq 35^{\circ}\text{C} \Delta \leq 0.2^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$
Water Temperature T1 (Control – Out) (3 checks)	$-0.5^{\circ}\text{C} \leq 6^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$ $-0.5^{\circ}\text{C} \leq 28^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$
Water Temperature T2 (Monitor – Out) (3 checks)	$-0.5^{\circ}\text{C} \leq 6^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$ $-0.5^{\circ}\text{C} \leq 28^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$ $-0.3^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.3^{\circ}\text{C}$
Water Temperature T3 (Monitor – In) (3 checks)	$-0.5^{\circ}\text{C} \leq 6^{\circ}\text{C} \Delta \leq 0.6^{\circ}\text{C}$ $-0.5^{\circ}\text{C} \leq 28^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$ $-0.5^{\circ}\text{C} \leq 40^{\circ}\text{C} \Delta \leq 0.5^{\circ}\text{C}$

Table 3.1. Verification Data

Section 4 – Maintenance

4.1 Recommended Maintenance

The Calibration and Test Unit contains no operator maintainable parts. Other than periodic cleaning of the exterior surfaces, no maintenance is required.

Cleaning instructions: Spray a mild detergent solution onto a towel and wipe down the unit. Take care to not saturate the towel with the solution, use only enough to dampen the unit. Towel dry the unit after wiping it down with the cleaning solution.

If a Mains fuse opens (inoperative heater), return the CTU for service. An open fuse is indicative of a problem that requires immediate service.

4.2 Calibration

The CTU should be returned to the manufacturer for calibration every two years. The date of last calibration is indicated by the international

symbol for calibration, an upside down solid black triangle, , which can be found on the top of the CTU.

4.3 Technical Assistance/Service

Customer Service: +1.844.823.5433, option 5, option 5,
8:30am - 6:00pm EST

Non-urgent Clinical and Technical Support: 877-267-2314

24-hour Clinical or Technical Support (during patient treatment): 866-840-9776

For assistance outside the US, contact your ARCTIC SUN™ Temperature Management System Distributor.

4.4 Disposal

Upon end of life, dispose of in accordance with local WEEE regulations or contact your local BARD® Supplier or Distributor and arrange for disposal.

APPENDIX A: Specifications

- Mains Input:
 - 115VAC, 60Hz, 4.5A
 - 230VAC, 50/60Hz, 2.25A
 - 100VAC, 50/60Hz, 5.0A
- Fuse: 6.3A Slow Blow (2)
- Power: 500W

APPENDIX B: Error Codes

Error Code	Description
Error 1	Pre-warm bypass flow error
Error 2	Pre-warm inlet pressure error
Error 3	Pre-warm nominal flow error
Error 4	Inlet pressure not stable
Error 5	Inlet pressure out of range
Error 11	Patient temperature 1 out of range at low test temperature (10°C)
Error 12	Patient temperature 1 out of range at mid-low test temperature (25°C)
Error 13	Patient temperature 1 out of range at mid-high test temperature (35°C)
Error 14	Patient temperature 1 out of range at high test temperature (40°C)
Error 15	Patient temperature 1 out of adjustment range limit
Error 16	Patient temperature 1 measurement unstable
Error 21	Patient temperature 2 out of range at low test temperature (10°C)
Error 22	Patient temperature 2 out of range at mid-low test temperature (25°C)
Error 23	Patient temperature 2 out of range at mid-high test temperature (35°C)
Error 24	Patient temperature 2 out of range at high test temperature (40°C)
Error 25	Patient temperature 2 out of adjustment range limit
Error 26	Patient temperature 2 measurement unstable
Error 30	Temperature out and patient temperature 2 difference out of range
Error 31	Temperature out measurement unstable
Error 32	Temperature out and simulator difference out of range
Error 40	Water temperature 1 is out of range at 6°C
Error 41	Water temperature 1 is out of range at 28°C
Error 42	Water temperature 1 is out of range at 40.5°C
Error 43	Water temperature 1 off of conversion chart table at 1°C
Error 44	Water temperature 1 off of conversion table at 47°C
Error 50	Water temperature 2 is out of range at 6°C
Error 51	Water temperature 2 is out of range at 28°C
Error 52	Water temperature 2 is out of range at 40.5°C
Error 53	Water temperature 2 off of conversion chart table at 1°C
Error 54	Water temperature 2 off of conversion table at 47°C
Error 60	Water temperature 3 is out of range at 6°C
Error 61	Water temperature 3 is out of range at 28°C
Error 62	Water temperature 3 is out of range at 40.5°C
Error 63	Water temperature 3 off of conversion chart table at 1°C
Error 64	Water temperature 3 off of conversion table at 47°C
Error 70	Water temperature 4 is out of range at 6°C
Error 71	Water temperature 4 is off of conversion table at 1°C
Error 80	Water check time out - Unable to reach calibration temperature
Error 81	Water check time out - Difference between water temperature and reference temperature is greater than 2°C
Error 82	Water check time out - other condition
Error 90	Heater test- flow is below 1 liter per minute
Error 91	Heater test- element 1 failure
Error 92	Heater test- element 2 failure
Error 93	Heater test- element 3 failure
Error 94	Heater test- element 4 failure
Error 95	Heater test- multiple elements fail
Error 99	Patient temperature out of calibration- no control
Error 100	Unable to initiate user diagnostic mode
Error 101	Unable to set calibration parameters
Error 102	Water level not full at start of calibration
Error 103	Switch settings incorrect
Error 104	Unknown error at startup
Error 105	Non recoverable navigation error
Error 106	Unknown error during pre-warm

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BD, the BD Logo, Bard, ArcticGel, Arctic Sun, Medivance and Simply Advanced are trademarks of Becton, Dickinson and Company or its affiliates. © 2025 BD. All rights reserved.
Medivance, Inc.
Becton Dickinson Ireland Ltd.