



BD Beaver™

Limbal Relaxing Incision Kit



Accurate Depth Knife – 550 or 600 micron

- Precise depth incisions with pre-set guard
- Bi-directional blade cuts left or right
- Multiple depth options available
- Single-use disposable knives



13 mm, Swivel Fixation Ring

- Guide marks on top of ring at 10° increments
- Less trauma than metal rings
- Single-use disposable swivel ring



Surgical Marking Pen

- For smooth, clear markings
- Prep resistant genietian violet ink
- Barrel legibly ruled in centimeters



Helping all people
live healthy lives

Advantages include:

- The precision of a BD Beaver™ Accurate Depth Knife
- The convenience of single-use disposable, high-quality products packaged together
- The elimination of stocking individual products for LRI incisions

Additional information:

Incision Nomogram* developed by Robert M. Kershner, MD
(www.seeadoc.com)

Incision Nomogram* developed by James P. Gills, MD
(www.stlukeseye.com)

Incision Nomogram* developed by Louis D. "Skip" Nichamin, MD
(www.eyenz.com)

*These nomograms are provided to you as objective sources of clinically relevant information for your reference. None of the material contained in this information was generated by BD.

BD Beaver™ Limbal Relaxing Incision Kits

Ordering Information:

| Catalog # | Description |
|-----------|--|
| 584055 | Limbal Relaxing Incision Kit (with 550 micron knife) |
| 584060 | Limbal Relaxing Incision Kit (with 600 micron knife) |
| 375055 | BD Beaver™ Accurate Depth Knife, 550 micron |
| 375060 | BD Beaver™ Accurate Depth Knife, 600 micron |
| 372242 | BD Beaver™ Swivel Fixation Ring |
| N/A | Surgical Marking Pen |

For ordering information call: **1-800-237-2174/1-201-847-7142**
or visit www.bd.com/ophthalmology

BD Medical
Ophthalmic Systems

1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/ophthalmology