

BD Beaver™ Family of Safety Knives for Ophthalmic Surgery



Helping all people
live healthy lives

Protection

Performance

Simplicity

The BD Beaver™ Family of Safety Knives for Ophthalmic Surgery

PROTECTION

The BD Beaver™ Safety Knives have an integrated retractable shield specifically designed to protect healthcare workers from accidental sharps injuries.

Safety Slit Knives

BD Beaver™ Xstar™ Safety Slit Knives

| Description | Single Bevel | Double Bevel |
|---------------------|--------------|--------------|
| 1.4 mm, angled 45° | 378214 | – |
| 1.5 mm, angled 45° | 378215 | – |
| 1.7 mm, angled 45° | 378217 | – |
| 2.2 mm, angled 45° | 378222 | – |
| 2.5 mm, angled 45° | 370145 | – |
| 2.6 mm, angled 45° | 370146 | – |
| 2.75 mm, angled 45° | 378227 | – |
| 2.8 mm, angled 45° | 378228 | 378828 |
| 3.0 mm, angled 45° | 378230 | 378830 |
| 3.2 mm, angled 45° | 378232 | 378832 |



BD Atomic Edge™ Safety Slit Knives

| Description | Knives |
|----------------------------------|--------|
| 3.0 mm, angled 45°, double bevel | 370157 |
| 2.8 mm, angled 45°, double bevel | 370159 |



Safety Sideport Knives

BD Beaver™ EdgeAhead™ Safety Sideport Knife

| Description | Knife |
|--------------|--------|
| .90 mm (20G) | 378231 |



BD Beaver™ Optimum™ Safety Sideport Knife

| Description | Knife |
|---------------|--------|
| 15°, straight | 378235 |



Safety Crescent Knives

BD Beaver™ Xstar™ Safety Crescent Knife

| Description | Knife |
|------------------------------|--------|
| 2.5 mm, angled 55°, bevel up | 378234 |



PERFORMANCE



The built-in shield protects the quality and consistency of the surgical blade from possible damage.

SIMPLICITY

The integrated safety shield is activated by the user through a spring-assisted slider that provides tactile and audible feedback in both the shielded and unshielded positions.


BD Beaver™ Safety Cataract Knife

Featuring Locking Retractable Shield




To Use

- Pull slider on handle to retract protective shield
- Retract shield to rear locked position to fully uncover blade



To Cover

- Push slider on handle to engage protective shield
- Advance shield to forward locked position at lock indicator line



To Dispose of

- Dispose of protected knife in an approved sharps container

For ordering information call:

1-866-513-SAFE/1-201-847-7142

or visit www.bd.com/safety

BD Beaver™ Family of Safety Knives for Ophthalmic Surgery

Protection

Performance

Simplicity



BD Medical
Ophthalmic Systems

1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/ophthalmology
1.866.513.SAFE (7233)
+1.781.906.7900