

BD Visitec™ LASEK



Helping all people
live healthy lives

*The convenience and innovation
you've come to expect from
BD Visitec, now for LASEK*



BD Visitec™ LASEK products continue with the convenience and innovation you have come to expect from the BD Visitec brand. Each single-use, sterile instrument utilizes a computer engineered steel component, matched to an ergonomically designed molded handle or hub.

Features and Benefits



585284

LASEK Dissector

Designed to create a uniform plane for segregating the epithelium from underlying tissue across the desired ablation area



585287

LASEK Cannula (Yee)

Olive tip designed to provide gentle manipulation with desired control, while multiple irrigation ports create a distribution of fluid when repositioning epithelium



585286

LASEK Microhoe (Yee)

Angled for optimum control with contoured face designed for precise resection of adherent epithelium



585285

LASEK Edge Lifter

Designed with a tapered tip and rounded sole to promote exact resection and prevent damage to underlying tissue when elevating the edge of the epithelial flap



BD, US Customer Service:

1 Becton Drive
Franklin Lakes, NJ 07417
tel: 800.237.2174 or 201.847.7142
fax: 800.804.7489 or 201.847.4833
special_services@BD.com
www.bd.com/ophthalmology

BD, UK Customer Service:

21 Between Towns Road
Cowley, Oxford OX4 3LY
England
tel: 44.1.865.781681
fax: 44.1.865.781627

Contact in other

European Countries:

11 rue Aristide Bergés-BP4
38800 Le Pont de Claix
France
tel: 33.4.76.68.32.27
fax: 33.4.76.68.32.87

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. ©2002 BD 587735-0203

Now combine your LASEK products with your favorite BD VISITEC™ products in a BD CUSTOMEYES™ surgical kit.

