


BD Atomic Edge™

Limbal Relaxing Incision Kit



Advantages include:

The precision of a BD Atomic Edge™ Accurate Depth Knife

The convenience of single-use high-quality products packaged together

The elimination of stocking individual products for LRI incisions

Kit Components

Accurate Depth Knife

- Equal in performance to diamond blades with single-use convenience
- Precise fixed blade exposure - Available in 550 and 600 micron
- Brand new sterile knife with each case



13 mm, Swivel Fixation Ring

- Guide marks on top of ring at 10° increments
- Less trauma than metal rings
- Single-use disposable swivel ring



Surgical Marking Pen

- For smooth, clear markings
- Prep resistant genietan violet risk
- Barrel legibly ruled in centimeters



BD Atomic Edge™ Limbal Relaxing Incision Kits

Cat. #	Description
586055	Limbal Relaxing Incision Kit (with 550 micron knife)
586060	Limbal Relaxing Incision Kit (with 600 micron knife)
370255	BD Atomic Edge™ Accurate Depth Knife, 550 micron
370260	BD Atomic Edge™ Accurate Depth Knife, 600 micron
372242	BD Beaver™ Swivel Fixation Ring
N/A	Surgical Marking Pen

Additional information:

Incision Nomogram* developed by James P. Gills, MD
(www.stlukeseye.com/professionals/lri_nomogram.htm)

*This nomogram is provided to you as an objective source of clinically relevant information for your reference. None of the material contained in this information was generated by BD.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2009 BD 587891-0409



Helping all people
live healthy lives

For ordering information call:
1.866.906.8080 option #2
or visit www.bd.com/safety

BD Medical
Ophthalmic Systems
1 Becton Drive
Franklin Lakes, NJ 07417
Tel: 800.237.2174 or
201.847.7142
www.bd.com/ophthalmology