

BD Visitec™ Single-Use Instruments

Meeting the Needs of Ophthalmic Surgery



Helping all people
live healthy lives

Protection Performance Simplicity

BD Medical - Ophthalmic Systems, the pioneer in developing single-use instruments for ophthalmic surgery, is proud to offer a wide variety of stainless steel instruments combining the precision you expect with the convenience of single-use.

Designed with safety in mind

Continuing a more than 70-year tradition of innovation, safety and unparalleled quality, we provide medical staff and patients high quality, innovative single-use instruments for ophthalmic surgery. Our single-use instruments are designed to reduce the risk of infection transmission and help you comply with today's safety requirements and recommendations.

BD Visitec™ Single-Use Instruments

Meeting the Needs of Ophthalmic Surgery

581409 Bonn Forceps (0.12 tip)

10 per box
100 mm long - 12 mm long tip.



581436 Moorfield Forceps

10 per box
100 mm long - 13 mm grooved section.



581422 Vannas Scissors

10 per box
120 mm.



581410 Capsulorhexis Forceps

10 per box
100 mm long - 12 mm long bent tip with iris stop.



581413 Eye Speculum

10 per box
40 mm x 30 mm - 14 mm open ears - 0.75 mm diameter steel.



581437 Westcott Scissors

10 per box
121 mm long - 22 mm scissor blade.



581411 Tying Forceps

10 per box
100 mm long - bent tip with micro-trays designed for 10/0 sutures.



581443 Eye Speculum

10 per box
40 mm x 30 mm - 14 mm open ears - 0.90 mm diameter steel.



581412 Buratto's Hook

10 per box
Angled 0.2 mm diameter, blunt tip.



581416 MacPherson Forceps

10 per box
100 mm long - 12 mm long tip.



581417 Needle Holder

10 per box
120 mm.



581423 Straight Micro-Manipulator

10 per box
Angled 0.35 mm diameter, blunt tip.



Protection ▶ A brand new, sterile instrument for each case ensures patient safety

Performance ▶ Precision stainless steel instruments with the convenience of single-use

Simplicity ▶ Easy to use, 100% traceability, eliminates the need for costly sterilization



BD Medical
Ophthalmic Systems
1 Becton Drive
Franklin Lakes, NJ 07417
Tel: 800.237.2174 or
201.847.7142
www.bd.com/ophthalmology