

REGULATED MEDICAL

Hazardous Waste

SHARPS

Infectious



PHARMACEUTICAL

Sharps
(No RCRA drugs)



WASTE COLLECTORS

Hazardous Waste

CHEMOTHERAPY

Trace



RCRA

P, U, D
listed drugs



SHARPS

With large access openings, including temporary and final closure, these sharps collectors are sturdy, stackable when full, and designed for easy handling and disposal.



PHARMACEUTICAL

These collectors feature sliding and dual opening hinged tops with VAT (Visual, Audible, Tactile) perimeter locks to collect and dispose unused, expired, or older pharmaceutical waste.



CHEMOTHERAPY

These collectors feature a wide entry port design for disposal of bulky chemotherapy waste and are designed to be leak resistant.



RCRA

These collectors are made from a minimum of 70% recycled material and are intended for solid RCRA hazardous waste.



BD XL Regulated Medical Waste Collectors Dimensions Guide

Reorder #	Product Size	Product Dimensions (H x W x D inches)
SHARPS		
305615	9 gal. Clear Hinge Top w/o gasket	18.75 x 17.75 x 11.75
305616	9 gal. Clear Slide Top w/o gasket	18.75 x 17.75 x 11.75
305601*	9 gal. Clear Hinge Top	18.75 x 17.75 x 11.75
305602*	9 gal. Clear Slide Top	18.75 x 17.75 x 11.75
305610	17 gal. Clear Hinge Top w/o gasket	23.5 x 20 x 14.75
305665*	17 gal. Clear Hinge Top	23.5 x 20 x 14.75
305609	19 gal. Clear Slide Top w/o gasket	26.25 x 20 x 14.75
305666*	19 gal. Clear Slide Top	26.25 x 20 x 14.75
PHARMACEUTICAL		
305634	9 gal. Hinge Top	18.75 x 17.75 x 11.75
305329	9 gal. Slide Top	18.75 x 17.75 x 11.75
305028	17 gal. Hinge Top	23.5 x 20 x 14.75
CHEMOTHERAPY		
305603	9 gal. Clear Hinge Top	18.75 x 17.75 x 11.75
305604	9 gal. Clear Slide Top	18.75 x 17.75 x 11.75
305614	17 gal. Clear Hinge Top	23.5 x 20 x 14.75
305613	19 gal. Clear Slide Top	26.25 x 20 x 14.75
RCRA		
305068	9 gal. Clear Hinge Top	18.75 x 17.75 x 11.75
305069	9 gal. Clear Slide Top	18.75 x 17.75 x 11.75
305071	17 gal. Clear Hinge Top	23.5 x 20 x 14.75
305072	19 gal. Clear Slide Top	26.25 x 20 x 14.75
TROLLEYS		
305556	9 gal. Standard Trolley	23.75 x 15.5 x 22.25
305567	9 gal. Foot Operated Trolley	26.5 x 15.5 x 22.25
305611	17 & 19 gal. Standard Trolley	26.5 x 17.5 x 23.75
305612	19 gal. Foot Operated Trolley	34.5 x 17.5 x 23.75

* Meets or exceeds US DOT packaging requirements 49CFR173.24, 49CFR173.24a, and 49CFR173.134.

For additional information
on BD products call
1-800-32-SHARP



BD Medical
4665 North Avenue
Oceanside, CA 92056
www.bd.com/sharps