

BD DISPOSAL SOLUTIONS ACCESSORIES



Helping all people
live healthy lives



Collector not included.

BD Recykleen™

Cabinets



BD Recykleen Cabinets are available in two sizes: small for use with 5.4 quart collectors and large cabinets for use with 2 and 3 gallon collectors.

Each size is available in both green and brown semi-translucent colors. Key product features include:

- An integrated window facilitating easy cleaning
- A universal lock providing security and preventing unwanted access to collectors
- The ability to clearly see the collector fill level
- A two-year warranty providing product replacement if it had manufacturer defects



Collectors not included.

Introducing BD Recykleen™ Cabinets and Trolleys

A new waste disposal management line made with recycled content.

BD Recykleen™

Trolleys



BD Recykleen Trolleys are available in standard and foot-operated versions and are used with large and extra large collectors.

Trolleys provide stability and mobility for large and extra large collectors, which can otherwise be very difficult to move when in use. BD Recykleen Trolleys also:

- Are designed to reduce the risk of collectors tipping over
- Enable the collector to be conveniently placed where it will be most used
- Come with locking wheels, ensuring the product is mobile only when you want it to be
- Provide the ability to open and close the collector for hands-free disposal (foot-operated trolleys only)
- Include a two-year warranty providing product replacement if it had manufacturer defects



Collector not included.



Collector not included.



Collector not included.



305099

305098

305096

305097

Collectors not included.

BD Recykleen™ Cabinets

Reference No.	Product Size	Qty/Case	Product Dimensions (H x W x D inches)
305096	5.4 Quart BD Recykleen Brown Cabinet	1	14 x 6 x 15
305097	2/3 Gallon BD Recykleen Brown Cabinet	1	15 x 7 x 19
305098	5.4 Quart BD Recykleen Green Cabinet	1	14 x 6 x 15
305099	2/3 Gallon BD Recykleen Green Cabinet	1	15 x 7 x 19



305091

305090

305092

305093

Collectors not included.

BD Recykleen™ Trolleys

Reference No.	Product Size	Qty/Case	Product Dimensions (H x W x D inches)
305090	9 Gallon Standard BD Recykleen Trolley	1	23.5 x 16.3 x 28.6
305092	17/19 Gallon Standard BD Recykleen Trolley	1	24.1 x 18 x 27.5
305091	9 Gallon Foot-operated BD Recykleen Trolley	1	23.5 x 16.3 x 28.6
305093	19 Gallon Foot-operated BD Recykleen Trolley	1	24.3 x 17.7 x 35.2

STABILIZERS



305021, 5/case
For all 5.4 Quart Collectors



305570, 5/case
For 8 Quart Multi-use Nestable Collectors



305552, 5/case
For 2/3 Gallon Counterbalanced Door Collectors

BRACKETS

305447, 10/case
Locking Bracket for all
5.4 Quart Collectors



305409, 10/case
Locking Bracket for 2/3 Gallon
Counterbalanced Door Collectors



305588, 6/case
Wire Bracket for 2/3 Gallon
Counterbalanced Door
Collectors



305485, 12/case
Bracket for Multi-use
Nestable, 8 & 14 Quart and
6 Gallon Collectors



305495, 12/case
XL Security Bracket for 3 Gallon
Plug Cap & 5 Gallon Chemo and
RCRA & 5 Gallon Multi-use One
Piece Collectors



305970, 12/case
Cup Bracket for Phlebotomy
Collectors



305497, 12/case
Quick Release Strap for 3 Gallon Plug Cap
& 5 Gallon and Multi-use One Piece
(except 1.4 qt) Collectors



305018, 6/case
Wire Bracket for 8 Quart
Multi-use Nestable Collectors



305496, 12/case
Security Bracket for
6.9 & 8.2 Quart Multi-use
One Piece Collectors



ACCESSORIES

OTHER



Collector not included.

**BD Recykleen™
Cabinets**



**BD RCRA
Collectors**



**BD Recykleen™
Large and XL Collectors**



**BD Recykleen™
Trolleys**

Contact your BD Sales Representative for an evaluation of your Sharps Collectors and Accessories needs and how BD Sharps can assist your department or institution.

You can also contact BD Sharps directly at **1-800-32-SHARP** for additional information on BD Disposal Solutions and Accessories or to request a visit from your BD Sales Representative.



BD Medical
4665 North Avenue
Oceanside, CA 92056
www.bd.com/sharps