


BD Eclipse™ Needle

Safety made simple

Easy to identify pink safety shield
fully encases needle upon activation

Wide, textured finger pad provides safe
activation area for thumb or forefinger

BD PrecisionGlide™ Needle Technology



Providing streamlined and effective injection safety

The BD Eclipse™ Needle is easy to use with minimal training, and supports streamlined and cost-effective option for every practitioner, every day*

- Allows for shielding of the needle tip immediately after injection, with a single finger stroke
- Activation requires minimal change in clinician technique for easy integration
- Specifically designed to reduce the risk of accidental needlestick injuries¹
- Audible click confirming activation to support your safety

Scan this code to watch a video on
how to use the **BD Eclipse™ Needle**:



* Within the BD safety hypodermic portfolio



BD Eclipse™ Needle

BD Eclipse™ Needles		
Product number	Description	Packaging
305757	30 G x ½ in.	100/box, 1,200/case
305758	27 G x ½ in.	100/box, 1,200/case
305767	25 G x 1 ½ in.	100/box, 1,200/case
305761	25 G x 1 in.	100/box, 1,200/case
305759	25 G x ⅝ in.	100/box, 1,200/case
305762	23 G x 1 in.	100/box, 1,200/case
303304	23 G x 1 ½ in.	100/box, 1,200/case
305765	21 G x 1 ½ in.	100/box, 1,200/case
305764	21 G x 1 in.	100/box, 1,200/case
305766	18 G x 1 ½ in.	100/box, 1,200/case

BD Eclipse™ Needles with detachable 1 mL BD Luer-Lok™ Syringe		
Product number	Description	Packaging
305778	30 G x ½ in.	50/box, 300/case
305789	27 G x ½ in.	50/box, 300/case
305780	25 G x ⅝ in.	50/box, 300/case

BD Eclipse™ Needles with detachable 3 mL BD Luer-Lok™ Syringe		
Product number	Description	Packaging
305787	25 G x 1 in.	50/box, 300/case
305781	25 G x ⅝ in.	50/box, 300/case
305782	23 G x 1 in.	50/box, 300/case
303306	23 G x 1 ½ in.	50/box, 300/case
305784	21 G x 1 ½ in.	50/box, 300/case
305779	21 G x 1 in.	50/box, 300/case

BD Eclipse™ Needle with detachable 5 mL BD Luer-Lok™ Syringe		
Product number	Description	Packaging
303307	23 G x 1 ½ in.	50/box, 300/case

BD Eclipse™ Needle with detachable 10 mL BD Luer-Lok™ Syringe		
Product number	Description	Packaging
303321	21 G x 1 ½ in.	50/box, 300/case

1. van der Molen F, Sluiter JK, Zwiderman KAH, Frings-Dresen MHW. Better effect of the use of a needle safety device in combination with an interactive workshop to prevent needle stick injuries. *Saf Sci.* 2011;49:1180-1186.

For more information on our products, please contact Medical Support Services at **800.555.7422** or medical.services@bd.com

To speak with your BD sales representative, call **888.237.2762**

bd.com

