

PRESCRIPTIVE INFORMATION

Use: For urological use only.

Caution

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- This product contains natural rubber latex which may cause allergic reactions.

Warnings

- Do not use ointments or lubricants having a petroleum base. This will damage latex. Visually inspect the product for any imperfections or surface deterioration prior to use.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Notes

- Single patient use, do not reuse.
- Sterile unless package is opened or damaged.